

Human Neutrophil Elastase Detection ELISA Kit

Catalog # 6057

For Research Use Only - Not Human or Therapeutic Use

PRODUCT SPECIFICATIONS

DESCRIPTION:	ELISA kit to quantify human neutrophil elastase
FORMAT:	Precoated 96-well ELISA Plate with removeable strips
ASSAY TYPE:	Sandwich ELISA
ASSAY TIME:	4 hours
STANDARD RANGE:	1000 - 16 pg/ml
NUMBER OF SAMPLES:	Up to 40 (duplicate) samples/plate
SAMPLE TYPES:	Liquid samples and biological fluids (pre-treatment acceptable)
RECOMMENDED SAMPLE DILUTIONS:	1:1 (at least)
CHROMOGEN:	TMB (read at 450 nm)
STORAGE:	-20°C for 12 months
VALIDATION DATA:	Intra-Assay (3.0-9.7%)/Inter-Assay (2.9-8.0%)/Spiking Test (94-95%)
NOTES:	N/A

Human Neutrophil Elastase Detection ELISA Kit

Catalog # 6057

For Research Use Only - Not Human or Therapeutic Use

INTRODUCTION

Neutrophils are frontline innate immune cells that defend the host from bacterial and viral infection through phagocytosis and the activity of intracellular antimicrobial proteins. Under certain pathological conditions or following *in-vitro* stimulation, neutrophils can release extracellular, web-like structures termed neutrophil extracellular traps (NETs). NETs are composed of chromatin decorated with DNA and a variety of neutrophil proteins, including neutrophil elastase (NE), myeloperoxidase (MPO), cathepsin G, and the S100A8/A9 (calprotectin) complex (1-3).

Neutrophil elastase is a serine protease stored in azurophilic (primary) granules and released upon neutrophil activation and degranulation. NE is a single-chain protein with a molecular weight of approximately 29 kDa and functions in extracellular matrix degradation, cytokine processing, and modulation of immune responses. Because of these activities, dysregulated NE contributes to the pathogenesis of acute lung injury, chronic obstructive pulmonary disease (COPD), cystic fibrosis, rheumatoid arthritis, sepsis, and ANCA-associated vasculitis (4,5).

For *in-vitro* studies of neutrophil activation and NET release, enzymatic activity assays (for example, NE and MPO activity) provide functional information about enzyme release and activation. Quantification of NE and MPO in cell culture medium by ELISA is widely used as a biochemical marker of NET release. Typical cell sources for NET assays include primary human neutrophils isolated from whole blood and established cell lines (e.g., HL-60 cells) (6).

Chondrex, Inc. provides a Human Neutrophil Elastase Detection ELISA Kit (Cat # 6057) suitable for assaying cell culture supernatant, biological fluid, serum, and plasma, as well as a Mouse MPO Detection ELISA Kit (Cat # 6051) that cross-reacts with human MPO. These products can be used to study NET release dynamics under neutrophil activation in experimental systems.

KIT COMPONENTS

Item	Quantity	Amount	Storage
Standard Human Neutrophil Elastase (60571)	1 vial	1000 pg, lyophilized	-20°C
Anti-Human Neutrophil Elastase Detection Antibody (60573)	1 vial	100 µl	-20°C
Solution C - Sample/Standard/Detection Antibody Dilution Buffer (30106)	1 bottle	50 ml	-20°C
Streptavidin Peroxidase (9029)	2 vials	50 µl	-20°C
Solution D - Streptavidin Peroxidase Dilution Buffer (9055)	1 bottle	20 ml	-20°C
TMB Solution (90023)	2 vials	0.2 ml	-20°C
Chromogen Dilution Buffer (90022)	1 bottle	20 ml	-20°C
Stop Solution - 2N Sulfuric Acid (9016)	1 bottle	10 ml	-20°C
Wash Buffer, 20X (9005)	1 bottle	50 ml	-20°C
Anti-Human Neutrophil Elastase Antibody Coated ELISA Plate (Brown)	1 each	96-well (8-well strips x 12)	-20°C

ASSAY OUTLINE




PLATE MAPPING

	1	2	3	4	5	6	7	8	9	10	11	12
A	1000	1000	S1	S1	S9	S9	S17	S17	S25	S25	S33	S33
B	500	500	S2	S2	S10	S10	S18	S18	S26	S26	S34	S34
C	250	250	S3	S3	S11	S11	S19	S19	S27	S27	S35	S35
D	125	125	S4	S4	S12	S12	S20	S20	S28	S28	S36	S36
E	63	63	S5	S5	S13	S13	S21	S21	S29	S29	S37	S37
F	31	31	S6	S6	S14	S14	S22	S22	S30	S30	S38	S38
G	16	16	S7	S7	S15	S15	S23	S23	S31	S31	S39	S39
H	B	B	S8	S8	S16	S16	S24	S24	S32	S32	S40	S40
	Standards		Samples									

NOTES BEFORE USING ASSAY

NOTE 1: It is recommended that the standard and samples be run in duplicate.

NOTE 2: Warm up all buffers to room temperature before use.

NOTE 3: Crystals may form in Wash Buffer, 20X when stored at cold temperatures. If crystals have formed, warm the wash buffer by placing the bottle in warm water until crystals are completely dissolved.

NOTE 4: Measure exact volume of buffers using a serological pipet, as extra buffer is provided.


NOTE 5: Cover the plate with plastic wrap or a plate sealer after each step to prevent evaporation from the outside wells of the plate.

NOTE 6: For partial reagent use, please see the assay protocol's corresponding step for the appropriate dilution ratio. For example, if the protocol dilutes 50 μ l of a stock solution in 10 ml of buffer for 12 strips, then for 6 strips, dilute 25 μ l of the stock solution in 5 ml of buffer. Partially used stock reagents may be kept in their original vials and stored at -20°C for use in a future assay.

NOTE 7: This kit contains animal components from non-infectious animals and should be treated as potential biohazards in use and for disposal.

ASSAY PROCEDURE

- Prepare Standard Dilutions:** The recommended standard range is 1000 - 16 pg/ml. Dissolve one vial of Standard (1000 pg/vial) in 1 ml of Sample/Standard Dilution Buffer (Solution C) and keep it as standard stock. Then serially dilute it with Solution C. For example, mix 250 μ l of the 1000 pg/ml solution with an equal volume of Solution C to make a 500 pg/ml solution, and then repeat it five more times for 250, 125, 63, 31, and 16 pg/ml standard solutions.



- Prepare Sample Dilutions:** Dilute samples at least 1:1 with Solution C depending on the estimated human neutrophil elastase levels in the samples. Two or three different sample dilutions are recommended if the human neutrophil elastase levels in the samples are unknown.
- Add Standards and Samples:** Add 100 μ l of standards, Solution C (blank), and samples to wells in duplicate. Incubate at room temperature for 2 hours.
- Wash:** Dilute 50 ml of 20X wash buffer in 950 ml of distilled water (1X wash buffer). Wash the plate with 1X wash buffer at least 3 times using a wash bottle with manifold or an automated plate washer. Empty the plate by inverting it and blotting on a paper towel to remove excess liquid. *Do not allow the plate to dry out.*

5. **Add Detection Antibody:** Dilute one vial of detection antibody in 10 ml Solution C. Add 100 μ l of detection antibody solution to each well and incubate at room temperature for 1 hour.

Strip #	Detection Antibody (μ l)	Solution C (ml)
2	17	1.7
4	33	3.3
6	50	5.0
8	66	6.6
10	82	8.2
12	100	10.0

6. **Wash:** Wash the plate with 1X wash buffer at least 3 times using a wash bottle with manifold or an automated plate washer. Empty the plate by inverting it and blotting on a paper towel to remove excess liquid. *Do not allow the plate to dry out.*
7. **Add Streptavidin Peroxidase Solution:** Prepare streptavidin peroxidase solution with Solution D as shown in the following table. Add 100 μ l of streptavidin peroxidase solution to each well and incubate at room temperature for 30 minutes.

Strip #	Streptavidin Peroxidase (μ l)	Solution D (ml)
2	8	1.7
4	17	3.3
6	25	5.0
8	33	6.6
10	42	8.2
12	50	10.0

8. **Wash:** Wash the plate with 1X wash buffer at least 3 times using a wash bottle with manifold or an automated plate washer. Empty the plate by inverting it and blotting on a paper towel to remove excess liquid. *Do not allow the plate to dry out.*
9. **Add TMB Solution:** Use new tubes when preparing TMB solution. Just prior to use, prepare TMB solution with Chromogen Dilution Buffer as shown in the following table. Add 100 μ l of TMB solution to each well immediately after washing the plate and incubate for 25 minutes at room temperature

Strip #	TMB (μ l)	Chromogen Dilution Buffer (ml)
2	34	1.7
4	66	3.3
6	100	5.0
8	132	6.6
10	164	8.2
12	200	10.0

10. **Stop:** Stop the reaction with 50 μ l of 2N Sulfuric Acid (Stop Solution) to each well.
11. **Read Plate:** Read the OD values at 450 nm. If the OD values of samples are greater than the OD values of the highest standard, re-assay the samples at a higher dilution. A 630 nm filter can be used as a reference.

CALCULATING RESULTS

1. Average the duplicate OD values for the standards, blanks (B), and test samples.
2. Subtract the averaged blank OD values from the averaged OD values of the standards and test samples.
3. Plot the OD values of the standards against the pg/ml of human neutrophil elastase standard. Using a log/log plot will linearize the data. Figure 1 shows an example of a standard curve where the standard range is 16 to 1000 pg/ml.
4. The pg/ml of human neutrophil elastase in test samples can be calculated using regression analysis. Multiply it by the sample dilution factor to obtain the human neutrophil elastase concentration (pg/ml) in the original test samples.




Figure 1 - A Typical Standard Curve for the Human Neutrophil Elastase Detection ELISA Kit

VALIDATION DATA

Table 1 - Reproducibility Data for the Human Neutrophil Elastase Detection ELISA Kit

Test	31 pg/ml	125 pg/ml	500 pg/ml
Intra-Assay CV (%)	9.7	6.7	3.0
Inter-Assay CV (%)	6.3	8.0	2.9
Spike Test* (%)	95%	94%	95%

*Known amounts of human neutrophil elastase were added to samples and then diluted with Sample/Standard Dilution Buffer (Solution C) to assay human neutrophil elastase by ELISA.

TROUBLESHOOTING

For frequently asked questions about assays and ELISAs, please see Chondrex, Inc.'s [ELISA FAQ](#) for more information.

REFERENCES

1. F. Castanheira, P. Kubes, Neutrophils and NETs in modulating acute and chronic inflammation. *Blood* **133**(20), 2178-2185 (2019).
2. B. van der Veen, M. P.J. de Winther, P. Heeringa, Myeloperoxidase: molecular mechanisms of action and their relevance to human health and disease. *Antioxid Redox Signal* **11**(11), 2899-937 (2009).
3. V. Papayannopoulos, K. Metzler, A. Hakkim, A. Zychlinsky, Neutrophil elastase and myeloperoxidase regulate the formation of neutrophil extracellular traps. *J Cell Biol* **191**(3), 677-91 (2010).
4. W. Lee, G. Downey, Leukocyte elastase: physiological functions and role in acute lung injury. *Am J Respir Crit Care Med* **164**(5), 896-904 (2001).
5. P. Henriksen, J. Sallénave, Human neutrophil elastase: mediator and therapeutic target in atherosclerosis. *Int J Biochem Cell Biol* **40**(6-7), 1095-100 (2008).
6. V. Papayannopoulos, K. D. Metzler, A. Hakkim, A. Zychlinsky, Neutrophil elastase and myeloperoxidase regulate the formation of neutrophil extracellular traps. *J. Cell Biol.* **191**, 677–691 (2010).