Immunization Grade Human Type || Collagen, Lyophilized

Catalog # 20051

For Research Use Only - Not Human or Therapeutic Use

DESCRIPTION:	Highly purified type II collagen (treated with 3M guanidine, DEAE-cellulose, and Na2HPO4)
APPLICATION:	Can be used for various purposes, such as a standard for collagen analysis by SDS-gel, substrate for MMP-1, MMP-8, and MMP-13, and as an antigen for immunizing animals for the induction of collagen-induced arthritis (CIA) in certain species and strains of experimental animals.
QUANTITY:	1 mg
FORM:	Lyophilized powder
SOURCE:	Human sternal cartilage
MOLECULAR WEIGHT:	300 kDa
PURITY:	>99% by SDS-PAGE, free of pepsin and proteoglycans
STORAGE:	4°C in the dark
STABILITY:	2 years
NOTES:	Soluble in acidic buffer (maximum 4 mg/ml), but difficult to dissolve in neutral buffer. In order to dissolve this product in a neutral buffer, first dissolve collagen in 0.01M or 0.05M acetic acid at 1-4 mg/ml and then add this collagen solution to the desired neutral buffer such as 2-10X Tris-NaCI (NaCI final concentration: 0.15-0.2M) or dialyze against 1X neutral buffer at 4°C.
	NOTE 1: To avoid fibril formation under neutral conditions, keep solution on ice.
	NOTE 2: The physicochemical property of type II collagen differs from type I collagen and type II collagen dissolved in a neutral buffer does not form a stable gel incubating at 37°C. In order to prepare a type II collagen gel, additional components such as type I collagen or other substances may be required.
REFERENCES:	R. Trelstad, et al. Biochemistry 9: 4993-4998 (1970)
	E. Miller, <i>Biochemistry</i> 10 : 1652-1659 (1971)
	D. Trentham, et al. J. Exp. Med. 146: 857-868 (1977)
	P. Wooley, et al. J. Immunology 135: 2443-2451 (1985)

© 2020 Chondrex, Inc. All Rights Reserved, 20051 2.0