

HMGB1 Detection ELISA kit

Catalog # 6010

For Research Use Only - Not Human or Therapeutic Use

PRODUCT SPECIFICATIONS

DESCRIPTION:	ELISA kit to quantify HMGB1
FORMAT:	96-well ELISA Plate with removeable strips
ASSAY TYPE:	Sandwich ELISA
ASSAY TIME:	2 hours*
STANDARD RANGE:	50 ng/ml to 0.8 ng/ml
NUMBER OF SAMPLES:	Up to 40 (duplicate) samples/plate
SAMPLE TYPES:	Cell culture, Serum, and Plasma
RECOMMENDED SAMPLE DILUTIONS:	1:1 (at least)
CHROMOGEN:	TMB (read at 450 nm)
STORAGE:	-20°C
VALIDATION DATA:	Human Serum: Intra-Assay (2.7-3.6%)/Inter-Assay (1.3-5.1%)/Recovery (101+/-25.1%) Mouse Serum: Intra-Assay (1.4-8.1%)/Inter-Assay (0.9-4.5%)/Recovery (104+/-8.3%)
NOTES:	*This kit has an overnight incubation step

HMGB1 Detection ELISA kit

Catalog # 6010

For Research Use Only - Not Human or Therapeutic Use

INTRODUCTION

HMGB1 (high mobility group box 1) (1) was recently rediscovered as a late lethal mediator of endotoxin (2) and is currently considered a pro-inflammatory cytokine that plays crucial roles in a variety of acute and chronic inflammatory diseases. HMGB1 contains 216 amino acids (3) and maintains 99% of its sequence identity among mice (4), rats (5), bovines (6), and humans (7). HMGB1 consists of three structural domains (8), termed “A box (9-85)”, “B box (88-162)”, and a negatively charged carboxyl terminus (186-216). Moreover, it has been previously shown that the B box recapitulates the pro-inflammatory activity whereas the A box acts as an antagonist of HMGB1 (9,10).


Several lines of evidence highlight the significance of HMGB1 in the immune inflammatory response. For example, it has been shown that HMGB1 is actively released by a variety of cells such as macrophages when stimulated by lipopolysaccharides (LPS), TNF- α , and IL-1 (2), and is passively released by injured or necrotic cells associated with collapsing cell structures. In fact, patients who died from septic shock had higher serum HMGB1 levels than surviving sepsis patients (11). Similarly, high serum HMGB1 levels are observed in sepsis animal models and in collagen-induced arthritis animal models (12). With regard to the function of the protein itself, HMGB1 has also been shown to stimulate the release of TNF- α and IL-1 (13,14), as well as bind LPS and synergistically increase peripheral blood mononuclear cell IL-6 production (15). Together, these observations demonstrate that HMGB1 plays important roles in the inflammatory cascade.

Chondrex, Inc. provides an HMGB1 Detection ELISA kit (Cat # 6010) to determine HMGB1 levels in cell culture media and sera. This kit contains enough reagents to measure 40 samples in duplicate together with standards.

KIT COMPONENTS


Item	Quantity	Amount	Storage
HMGB1 Standard (60101)	2 vials	50 μ l/vial	-20°C
Capture Antibody (Anti-HMGB1 Monoclonal Antibody) (60102)	1 vial	100 μ l/vial	-20°C
Detection Antibody (Anti-HMGB1 Monoclonal Antibody) (60103)	1 vial	Lyophilized	-20°C
Solution A - Coating Buffer (9052)	1 bottle	10 ml	-20°C
Solution B - Sample/Standard Dilution Buffer (601010)	1 bottle	20 ml	-20°C
Solution C - Detection Antibody Dilution Buffer (60106)	1 bottle	10 ml	-20°C
Solution D - Streptavidin Peroxidase Dilution Buffer (9055)	1 bottle	20 ml	-20°C
Streptavidin Peroxidase (9029)	2 vials	50 μ l	-20°C
TMB Solution (contains DMSO) (90023)	2 vials	0.2 ml	-20°C
Chromogen Dilution Buffer (90022)	1 bottle	20 ml	-20°C
Stop Solution - 2N Sulfuric Acid (9016)	1 bottle	10 ml	-20°C
Wash Buffer, 20X (9005)	1 bottle	50 ml	-20°C
ELISA Plate	1 each	96-well (8-well strips x 12)	-20°C

ASSAY OUTLINE



* Use one tip per sample or standards. Do not cross-contaminate samples or standards by re-using pipet tips. A multi-channel pipet is recommended.

PLATE LAYOUT



NOTES BEFORE USING ASSAY

NOTE 1: It is recommended that the standard and samples be run in duplicate.

NOTE 2: Warm up all buffers to room temperature before use.

NOTE 3: Crystals may form in Wash Buffer, 20X when stored at cold temperatures. If crystals have formed, warm the wash buffer by placing the bottle in warm water until crystals are completely dissolved.

NOTE 4: Measure exact volume of buffers using a serological pipet, as extra buffer is provided.

NOTE 5: Cover the plate with plastic wrap or a plate sealer after each step to prevent evaporation from the outside wells of the plate.

NOTE 6: For partial reagent use, please see the assay protocol's corresponding step for the appropriate dilution ratio. For example, if the protocol dilutes 50 μ l of a stock solution in 10 ml of buffer for 12 strips, then for 6 strips, dilute 25 μ l of the stock solution in 5 ml of buffer. Partially used stock reagents may be kept in their original vials and stored at -20°C for use in a future assay.

NOTE 7: This kit contains animal components from non-infectious animals and should be treated as potential biohazards in use and for disposal.


NOTE 8: This kit can be used to determine HMGB1 levels in sera and cell culture media samples. However, special concern should be considered for assaying HMGB1 in human serum because autoantibodies to HMGB1 are determined in 9-89% of sera from patients with autoimmune and inflammatory diseases (16-19). These reports indicate that human serum polyclonal antibodies to HMGB1 might mask the epitopes recognized by the capture and detection antibodies used in this kit, resulting in interference against the assay.

ASSAY PROCEDURE

1. **Add Capture Antibody:** Dilute 1 vial (100 μ l) of Capture Antibody with 10 ml of Coating Buffer (Solution A). Add 100 μ l of capture antibody solution to each well and incubate at 4°C overnight. If planning to use less, the remaining stock solution can be stored in its original vial at -20°C.

Strip #	Capture Antibody (μ l)	Solution A (ml)
2	17	1.7
4	33	3.3
6	50	5.0
8	66	6.6
10	82	8.2
12	100	10.0

2. **Prepare Standard Dilutions:** The recommended standard range is 0.8-50 ng/ml. Dilute one vial of HMGB1 Standard with 950 μ l of Sample/Standard Dilution Buffer (Solution B) - 100 ng/ml. Prepare serial dilutions of the standard by mixing 200 μ l of the 100 ng/ml standard with 200 μ l of Solution B - 50 ng/ml. Then repeat this procedure to make six more serial dilutions of standard - 25, 12.5, 6.25, 3.1, 1.6, and 0.8 ng/ml solutions. Partially used 100 ng/ml standard stock **cannot** be saved for future assays. Discard unused, diluted standard solution. Chondrex, Inc. recommends making fresh standard and serial dilutions for each assay.



- Prepare Sample Dilutions:** Centrifuge samples at 10,000 rpm at 4°C for 3 minutes to remove insoluble materials and lipids and use the supernatant as samples. Dilute samples at least 1:1 with an equal volume of Solution B. For example, take 100 µl of a serum, and mix with 100 µl of Solution B.
- Prepare Detection Antibody:** Reconstitute one vial of Detection Antibody with 50 µl of distilled water to make a stock solution. Dilute the 50 µl of detection antibody stock solution in 5 ml of Detection Antibody Dilution Buffer (Solution C). If planning to use less, the remaining stock solution can be stored in its original vial at -20°C.

Strip #	Detection Antibody (µl)	Solution C (ml)
2	8	0.8
4	17	1.7
6	25	2.5
8	33	3.3
10	42	4.2
12	50	5.0

- Dilute Wash Buffer:** Dilute 50 ml of 20X wash buffer in 950 ml of distilled water (1X wash buffer). Wash the plate with 1X wash buffer at least 3 times using a wash bottle with manifold or an automated plate washer. Empty the plate by inverting it and blotting on a paper towel to remove excess liquid. *Do not allow the plate to dry out.*
- Add Standards, Samples, and Detection Antibody:** Vortex standards, samples, and detection antibody tubes well. Add 50 µl of Solution B (blank), standards, and samples to appropriate wells. Add 50 µl of diluted detection antibody solution to all wells. Mix all wells by pipetting or use a plate shaker. Cover the plate with a plate sealer and incubate at 37°C for 1 hour, then transfer plate to 4°C overnight.
- Wash:** Wash the plate with 1X wash buffer at least 3 times using a wash bottle with manifold or an automated plate washer. Empty the plate by inverting it and blotting on a paper towel to remove excess liquid. *Do not allow the plate to dry out.*
- Add Streptavidin Peroxidase:** Dilute one vial of Streptavidin Peroxidase in 10 ml of Streptavidin Peroxidase Dilution Buffer (Solution D). If planning to use less, the remaining stock solution can be stored in its original vial at -20°C. Add 100 µl of streptavidin peroxidase solution to each well and incubate at room temperature for 30 minutes.

Strip #	Streptavidin Peroxidase (µl)	Solution D (ml)
2	9	1.8
4	17	3.4
6	25	5.0
8	33	6.6
10	41	8.2
12	50	10.0

9. **Wash:** Wash the plate with 1X wash buffer at least 3 times using a wash bottle with manifold or an automated plate washer. Empty the plate by inverting it and blotting on a paper towel to remove excess liquid. *Do not allow the plate to dry out.*
10. **Add TMB:** Use new tubes when preparing TMB. Dilute one vial of TMB with 10 ml of Chromogen Dilution Buffer just prior to use. Add 100 μ l of TMB solution to all wells immediately after washing the plate. Incubate for 30 minutes at room temperature. If planning to use less, the remaining stock solution can be stored in its original vial at -20°C .


Strip #	TMB (μ l)	Chromogen Dilution Buffer (ml)
2	34	1.7
4	66	3.3
6	100	5.0
8	132	6.6
10	164	8.2
12	200	10.0

11. **Stop:** Add 50 μ l of 2N sulfuric acid (Stop Solution) to each well.
12. **Read Plate:** Read the OD values at 450 nm (a 630 nm filter can be used as a reference) immediately. If the OD values of samples are greater than the OD values of the highest standard, re-assay the samples at a higher dilution.

CALCULATING RESULTS

1. Average the duplicate OD values for the blank, standards, and samples.
2. Subtract the averaged blank (B) OD value from the averaged standard and sample OD values.
3. Plot the OD values of standards against the concentration of HMGB1 (ng/ml) using a log scale. Figure 1 shows a typical standard curve where the HMGB1 range is from 0.8-50 ng/ml.
4. The concentration of HMGB1 (ng/ml) in samples can be calculated using regression analysis. Multiply the results by the dilution factors (usually 2 without extra dilution). For additional assistance, please download a [sample calculation worksheet](#) from www.chondrex.com.

Figure 2 - A Typical Standard Curve for the HMGB1 Detection ELISA Kit



VALIDATION DATA

Table 1 - Reproducibility Data for the HMGB1 Detection ELISA Kit

Human Serum

Test At	2 ng/ml	7.5 ng/ml	30 ng/ml
Intra-Assay CV (%)	2.7	3.4	3.6
Inter-Assay CV (%)	1.3	5.1	3.1

Mouse Serum

Test At	2 ng/ml	7.5 ng/ml	30 ng/ml
Intra-Assay CV (%)	1.4	3.2	8.1
Inter-Assay CV (%)	3.7	4.5	0.9

Recovery

Species	Averaged Recovery Results
Human (5 sera)	101 ± 25.1%
Mouse (5 sera)	104 ± 8.3%

Specificity: Average cross reactivity with bovine HMGB2 is 12.4%.

TROUBLESHOOTING

For frequently asked questions about assays and ELISAs, please see Chondrex, Inc.'s [ELISA FAQ](#) for more information.

REFERENCES

1. L. Einck, M. Bustin, The Intracellular Distribution and Function of the High Mobility Group Chromosomal Proteins. *Exp Cell Res* **156**, 295-310 (1985).
2. H. Wang, O. Bloom, M. Zhang, J. Vishnubhakat, M. Ombrellino, *et al.*, HMG-1 as a Late Mediator of Endotoxin Lethality in Mice. *Science* **285**, 248-51 (1999).
3. J. Parkkinen, E. Raulo, J. Merenmies, R. Nolo, E. Kajander, *et al.*, Amphoterin, the 30-kDa Protein in a Family of HMG1-type Polypeptides. Enhanced Expression in Transformed Cells, Leading Edge Localization, and Interactions With Plasminogen Activation. *J Biol Chem* **268**, 19726-38 (1993).
4. S. Ferrari, L. Ronfani, S. Calogero, M. Bianchi, The Mouse Gene Coding for High Mobility Group 1 Protein (HMG1). *J Biol Chem* **269**, 28803-8 (1994).
5. G. Paonessa, R. Frank, R. Cortese, Nucleotide Sequence of Rat Liver HMG1 cDNA. *Nucleic Acids Res* **15**, 9077 (1987).
6. D. Kaplan, C. Duncan, Full Length cDNA Sequence for Bovine High Mobility Group 1 (HMG1) Protein. *Nucleic Acids Res* **16**, 10375 (1988).

7. L. Wen, J. Huang, B. Johnson, G. Reeck, A Human Placental cDNA Clone That Encodes Nonhistone Chromosomal Protein HMG-1. *Nucleic Acids Res* **17**, 1197-214 (1989).
8. J. Li, H. Wang, J. Mason, J. Levine, M. Yu, *et al.*, Recombinant HMGB1 With Cytokine-Stimulating Activity. *J Immunol Methods* **289**, 211-23 (2004).
9. M. Bustin, Regulation of DNA-dependent Activities by the Functional Motifs of the High-Mobility-Group Chromosomal Proteins. *Mol Cell Biol* **19**, 5237-46 (1999).
10. H. Yang, M. Ochani, J. Li, X. Qiang, M. Tanovic, H. Harris, *et al.* Reversing established sepsis with antagonists of endogenous high-mobility group box 1. *Proc Natl Acad Sci U S A.* **101**:296-301 (2004).
11. U. Andersson, H. Wang, K. Palmblad, A. Aveberger, O. Bloom, H. Erlandsson-Harris, *et al.* High mobility group 1 protein (Hmg-1) stimulates proinflammatory cytokine synthesis in human monocytes. *J Exp Med.* **192**:565-570 (2000).
12. R. Kokkola, J. Li, E. Sundberg, A. Aveberger, K. Palmblad, *et al.*, Successful Treatment of Collagen-Induced Arthritis in Mice and Rats by Targeting Extracellular High Mobility Group Box Chromosomal Protein 1 Activity. *Arthritis Rheum* **48**, 2052-8 (2003).
13. H. Wang, J. Vishnubakat, O. Bloom, M. Zhang, M. Ombrellino, *et al.*, Proinflammatory Cytokines (Tumor Necrosis Factor and Interleukin 1) Stimulate Release of High Mobility Group protein-1 by Pituicytes. *Surgery* **126**, 389-92 (1999).
14. E. Abraham, J. Arcaroli, A. Carmody, H. Wang, K. Tracey, HMG-1 as a Mediator of Acute Lung Inflammation. *J Immunol* **165**, 2950-4 (2000).
15. H. Hreggvidsdottir, T. Ostberg, H. Wähämaa, H. Schierbeck, A. Aveberger, *et al.*, The Alarmin HMGB1 Acts in Synergy With Endogenous and Exogenous Danger Signals to Promote Inflammation. *J Leukoc Biol* **86**, 655-62 (2009).
16. B. Wittemann, G. Neuer, H. Michels, H. Truckenbrodt, F. Bautz, Autoantibodies to Nonhistone Chromosomal Proteins HMG-1 and HMG-2 in Sera of Patients With Juvenile Rheumatoid Arthritis. *Arthritis Rheum* **33**, 1378-83 (1990).
17. H. Uesugi, S. Ozaki, J. Sobajima, F. Osakada, H. Shirakawa, *et al.*, Prevalence and Characterization of Novel pANCA, Antibodies to the High Mobility Group Non-Histone Chromosomal Proteins HMG1 and HMG2, in Systemic Rheumatic Diseases. *J Rheumatol* **25**, 703-9 (1998).
18. J. Sobajima, S. Ozaki, H. Uesugi, F. Osakada, M. Inoue, *et al.*, High Mobility Group (HMG) Non-Histone Chromosomal Proteins HMG1 and HMG2 Are Significant Target Antigens of Perinuclear Anti-Neutrophil Cytoplasmic Antibodies in Autoimmune Hepatitis. *Gut* **44**, 867-73 (1999).
19. J. Sobajima, S. Ozaki, F. Osakada, H. Uesugi, H. Shirakawa, *et al.*, Novel Autoantigens of Perinuclear Anti-Neutrophil Cytoplasmic Antibodies (P-ANCA) in Ulcerative Colitis: Non-Histone Chromosomal Proteins, HMG1 and HMG2. *Clin Exp Immunol* **107**, 135-40 (1997).