

DNA Assay Kit

Catalog # 6023

For Research Use Only - Not Human or Therapeutic Use

PRODUCT SPECIFICATIONS

DESCRIPTION:	Assay kit to quantify DNA
FORMAT:	96-well ELISA plate with non-removeable strips
ASSAY TYPE:	Dye-based assay/Fluorescence-based assay
ASSAY TIME:	<0.5 hour
STANDARD RANGE:	10 µg/ml to 0.16 µg/ml
NUMBER OF SAMPLES:	Up to 40 (duplicate) samples/plate
SAMPLE TYPES:	Solubilized tissue & cell culture samples (pre-treatment acceptable)
RECOMMENDED SAMPLE DILUTIONS:	Varies (Typical DNA levels are usually 0.3-5 µg/mg wet tissue and 1-5 µg/1x10 ⁶ cells)
FLUORESCENCE:	Excitation 360 nm/Emission 460 nm
STORAGE:	-20°C for 12 months
VALIDATION DATA:	Intra-Assay (1.7-7.6%)/Inter-Assay (5.2-8.0%)/Spiking Test (105.1-116.2%)
NOTES:	

DNA Assay Kit

Catalog # 6023

For Research Use Only - Not Human or Therapeutic Use

INTRODUCTION

Deoxyribonucleic acid (DNA) is a distinct cell component, therefore DNA amount is correlated with cell number (1). For example, in cartilage tissue engineering, artificial cartilage quality is evaluated by DNA amounts translated as chondrocyte numbers, as well as amounts of collagen and glycosaminoglycans (GAGs) in extracellular matrix components (ECM) (2, 3). DNA can be quantified by the ratio of absorbance at 260 nm and 280 nm, but this method is affected by contaminating proteins, RNA, and chemicals in the sample, requiring additional DNA isolation steps. For tissue analysis, Chondrex, Inc. provides an alternative DNA assay kit employing the Hoechst 33258 fluorescent dye which specifically binds to Adenine-Thymine base pairs, resulting in fluorescence at excitation 360 nm/emission 460 nm. Because the dye-DNA binding and the fluorescence intensity are unaffected by contaminating proteins and other substances in an optimized assay condition, this DNA assay kit works accurately with samples containing other analytes.

In order to facilitate further studies in this field, Chondrex, Inc. also provides various native collagen detection ELISAs, a Hydroxyproline Assay Kit (Cat # 6017), Sirius Red/ Fast Green Staining Kit (Cat # 9046), Sirius Red Total Collagen Detection Kit (Cat # 9062), and GAGs assay kit (Cat # 6022). Please visit www.chondrex.com or contact support@chondrex.com for more information.

KIT COMPONENTS

Item	Quantity	Amount	Storage
Calf Thymus DNA Standard (60231)	1 vial	100 µg/ml, 100 µl	-20°C
Reaction Solution (60232)	1 bottle	50 ml	-20°C
Detection Solution Concentrate 200X (60233)	1 vial	50 µl	-20°C
Black ELISA Plate	1 each	96-well	-20°C

ASSAY OUTLINE





PLATE MAPPING



NOTES BEFORE USING ASSAY

NOTE 1: It is recommended that the standard and samples be run in duplicate.

NOTE 2: Warm up all buffers to room temperature before use.

NOTE 3: Crystals may form in Wash Buffer, 20X when stored at cold temperatures. If crystals have formed, warm the wash buffer by placing the bottle in warm water until crystals are completely dissolved.

NOTE 4: Measure exact volume of buffers using a serological pipet, as extra buffer is provided.


NOTE 5: Cover the plate with plastic wrap or a plate sealer after each step to prevent evaporation from the outside wells of the plate.

Table 1 - Solubilizing Reagents in the DNA Assay Kit

Reagents	Concentration in Samples
Guanidine	0.5 M
SDS	0.025%
Tween 20	0.0025%
Triton X - 100	0.025%

ASSAY PROCEDURE

- Prepare Standard Dilutions:** Take 40 μl of DNA Standard Solution and mix with 360 μl of Reaction Solution (10 $\mu\text{g}/\text{ml}$). Then serially dilute it with Reaction Solution. For example, mix 200 μl of the standard (10 $\mu\text{g}/\text{ml}$) with an equal volume of Reaction Solution to make a 5 $\mu\text{g}/\text{ml}$ solution, and then repeat it five more times for 2.5, 1.3, 0.63 0.31, and 0.16 $\mu\text{g}/\text{ml}$ solutions. Chondrex, Inc. recommends making fresh serial dilutions for each assay.



- Prepare Sample Dilutions:** Dilute samples with Reaction Solution. The optimal sample dilution depends on the sample and preparation methods; therefore, Chondrex, Inc. recommends multiple dilutions as a preliminary study.
- Prepare 1X Detection Solution:** Mix 5 μl of 200X Concentrated Detection Solution and 1 ml Reaction Solution for each well. For example, 8 samples, 7-point standards, one blank (all in duplicate) will require 1.6 ml of the 1X Detection Solution. Mix 8 μl of 200X Concentrated Detection Solution and 1.6 ml Reaction Solution.

NOTE: Prepare the 1X solution just before use. Do not store and reuse the mixed solution for the next assay.

- Add Standards and Samples:** Add 50 μl of diluted standards and samples into wells.
- Add 1X Detection Solution:** Add 50 μl of 1X detection solution into all wells, and then incubate at room temperature for 5 minutes.
- Read plate:** Read the plate at excitation 360 nm/emission 460 nm.


CALCULATING RESULTS

- Average the duplicate fluorescence intensity (FI) values for the blank, standards, and test samples.
- Subtract the averaged "blank" (B) FI values from the averaged FI values of standards and test samples.
- Plot the FI values of the standards against the DNA concentration ($\mu\text{g}/\text{ml}$). Using a log/log plot will linearize the data. Figure 1 shows an example of a DNA standard curve from 0.16-10 $\mu\text{g}/\text{ml}$.
- The $\mu\text{g}/\text{ml}$ of DNA in test samples can be calculated using regression analysis. Multiply it by the sample dilution factor to obtain the DNA concentration ($\mu\text{g}/\text{ml}$) in the original sample.

NOTE 1: DNA levels in samples depend on the sample solubilization efficiency of individual protocols consequently, converting DNA levels to cell numbers is not always applicable. Therefore, DNA levels as correlating to cell numbers can be used to compare samples including positive and negative controls.

NOTE 2: To correlate DNA amounts to cell numbers, validating your DNA solubilization protocol is required. Please refer to the references for more information.

Figure 1 - A Typical Standard Curve for the DNA Assay



ASSAY VALIDATION

Table 2 - Reproducibility Data for the DNA Assay Kit

Test	0.31 µg/ml	0.13 µg/ml	5 µg/ml
Intra-Assay CV (%)	7.6	3.6	1.7
Inter-Assay CV (%)	7.3	5.2	8.0
Spike Test* (%)	115.1%	116.2%	105.1%

* Known amounts of DNA were added to samples and then diluted with Reaction Solution for assaying DNA.

TROUBLESHOOTING

For frequently asked questions about assays and ELISAs, please see Chondrex, Inc.'s [Assay FAQ](#) for more information.

REFERENCES

1. V. M. Quent, D. Loessner, T. Friis, J. C. Reichert, D. W. Huttmacher, Discrepancies between metabolic activity and DNA content as tool to assess cell proliferation in cancer research. *J Cell Mol Med* **14**, 1003-1013 (2010).
2. J. H. Yoon, J. Halper, Tendon proteoglycans: biochemistry and function. *J Musculoskelet Neuronal Interact* **5**, 22-34 (2005).
3. J. A. Buckwalter, H. J. Mankin, Articular cartilage: tissue design and chondrocyte-matrix interactions. *Instr Course Lect* **47**, 477-486 (1998).