

# Mouse Mast Cell Protease-1 (MCPT-1) Detection ELISA Kit

Catalog # 6046

*For Research Use Only - Not Human or Therapeutic Use*

## PRODUCT SPECIFICATIONS

|                               |  |
|-------------------------------|--|
| DESCRIPTION:                  | ELISA kit to quantify mouse mast cell protease-1 (MCPT-1)            |
| FORMAT:                       | Precoated 96-well ELISA Plate with removeable strips                 |
| ASSAY TYPE:                   | Sandwich ELISA   |
| ASSAY TIME:                   | 4 hours  |
| STANDARD RANGE:               | 0.8 - 50 ng/ml   |
| NUMBER OF SAMPLES:            | Up to 40 (duplicate) samples/plate                                   |
| SAMPLE TYPES:                 | Mouse serum, culture medium, and biological fluids                   |
| RECOMMENDED SAMPLE DILUTIONS: | 1:1 (at least)   |
| CHROMOGEN:                    | TMB (read at 450 nm)   |
| STORAGE:                      | -20°C for 12 months  |
| VALIDATION DATA:              | Intra-Assay (3.1-6.7%)/Inter-Assay (1.7-8.0%)/Spiking Test (91-108%) |
| NOTES:                        | N/A  |

# Mouse Mast Cell Protease-1 (MCPT-1) Detection ELISA Kit

Catalog # 6046

*For Research Use Only - Not Human or Therapeutic Use*

## INTRODUCTION

Mast cells (MCs) are primarily recognized for their roles in allergic diseases, although they also participate in other health and disease processes. These cells accumulate at the sites of Th2 cell activation and initiate immediate hypersensitivity reactions. During IgE-associated biological responses, antigen-specific IgE antibodies bound to FcεRI on the plasma membrane of MCs undergo antigen-dependent cross-linking, leading to the secretion of biologically active products involved in allergic reactions. These products include vasoactive amines (histamines), neutral proteases (chymases), proteoglycans (heparans), and various cytokines and growth factors through a process called degranulation(1).

Mouse mast cell protease-1 (MCPT-1) is a β-chymase, a type of serine protease stored and secreted by intestinal mucosal mast cells, which are found in the intestinal epithelium (2). It shares 74% of its amino acid homology with its rat counterpart, rat mast cell protease-II (MCP-II) (3).


MCPT-1 serves as a marker for mast cell activation and degranulation. Although it is expressed constitutively and is detectable in the sera of normal mice, parasites in the gut cause systemic levels to increase dramatically within two days and peak at two weeks following infection (4). Deficiency in MCPT-1 is associated with significantly delayed expulsion of *Trichinella spiralis*, indicating its important role in the host defense against intestinal parasites (5). Elevated mouse MCPT-1 levels are also observed during intestinal allergic hypersensitivity reactions and have been reported as a marker for food allergy models (6). The mechanism is not fully understood, although it has been shown to be a tight junction-breaking protease, thus increasing the permeability of the small intestine (7).

Chondrex, Inc. provides a mouse MCPT-1 Detection ELISA kit (Cat # 6046) to determine MCP-1 levels in biological fluids and sera. This kit contains enough reagents to measure 40 samples in duplicate together with standards. Please visit [www.chondrex.com](http://www.chondrex.com) or contact [support@chondrex.com](mailto:support@chondrex.com) for more information.


## KIT COMPONENTS

| Item  | Quantity | Amount                       | Storage |
|---|----------|------------------------------|---------|
| Standard mouse MCPT-1 (60461)   | 1 vial   | 50 ng, lyophilized           | -20°C   |
| Anti-mouse MCPT-1 Detection Antibody (60463)                            | 1 vial   | 100 µl                       | -20°C   |
| Solution B - Sample/Standard/Detection Antibody Dilution Buffer (67015) | 1 bottle | 50 ml                        | -20°C   |
| Solution D - Streptavidin Peroxidase Dilution Buffer (9055)             | 1 bottle | 20 ml                        | -20°C   |
| Streptavidin Peroxidase (9029)  | 2 vials  | 50 µl                        | -20°C   |
| TMB Solution (90023)  | 2 vials  | 0.2 ml                       | -20°C   |
| Chromogen Dilution Buffer (90022)                                       | 1 bottle | 20 ml                        | -20°C   |
| Stop Solution - 2N Sulfuric Acid (9016)                                 | 1 bottle | 10 ml                        | -20°C   |
| Wash Buffer, 20X (9005)   | 1 bottle | 50 ml                        | -20°C   |
| Anti-mouse MCP-1 Antibody Coated ELISA Plate (Red)                      | 1 each   | 96-well (8-well strips x 12) | -20°C   |

## ASSAY OUTLINE



## PLATE MAPPING



## NOTES BEFORE USING ASSAY

NOTE 1: It is recommended that the standard and samples be run in duplicate.

NOTE 2: Warm up all buffers to room temperature before use.

NOTE 3: Crystals may form in Wash Buffer, 20X when stored at cold temperatures. If crystals have formed, warm the wash buffer by placing the bottle in warm water until crystals are completely dissolved.

NOTE 4: Measure exact volume of buffers using a serological pipet, as extra buffer is provided.


NOTE 5: Cover the plate with plastic wrap or a plate sealer after each step to prevent evaporation from the outside wells of the plate.

NOTE 6: For partial reagent use, please see the assay protocol's corresponding step for the appropriate dilution ratio. For example, if the protocol dilutes 50  $\mu$ l of a stock solution in 10 ml of buffer for 12 strips, then for 6 strips, dilute 25  $\mu$ l of the stock solution in 5 ml of buffer. Partially used stock reagents may be kept in their original vials and stored at -20°C for use in a future assay.

NOTE 7: This kit contains animal components from non-infectious animals and should be treated as potential biohazards in use and for disposal.

## ASSAY PROCEDURE

- Prepare Standard Dilutions:** The recommended standard range is 0.8 - 50 ng/ml. Dissolve one vial of Standard (50 ng/vial) in 1 ml of Sample/Standard/Detection Antibody Dilution Buffer (Solution B) and keep it as a standard stock. Then serially dilute it with Solution B. For example, mix 250  $\mu$ l of the 50 ng/ml solution with an equal volume of Solution B to make a 25 ng/ml solution, and then repeat it five more times for 12.5, 6.3, 3.1, 1.6, and 0.8 ng/ml standard solutions.



- Prepare Sample Dilutions:** Dilute samples at least 1:1 with Solution B depending on the estimated mMCP-1 level in the samples. Two or three different sample dilutions are recommended if the mMCP-1 levels in the samples are unknown.
- Add Standards and Samples:** Add 100  $\mu$ l of standards, Solution B (blank), and samples to wells in duplicate. Incubate at room temperature for 2 hours.
- Wash:** Dilute 50 ml of 20X wash buffer in 950 ml of distilled water (1X wash buffer). Wash the plate with 1X wash buffer at least 3 times using a wash bottle with manifold or an automated plate washer. Empty the plate by inverting it and blotting on a paper towel to remove excess liquid. *Do not allow the plate to dry out.*

5. **Add Detection Antibody:** Dilute one vial of Detection Antibody in 10 ml Sample/Standard/Antibody Dilution buffer (Solution B). Add 100  $\mu$ l of detection antibody solution to each well and incubate at room temperature for 1 hour.

| Strip # | Detection Antibody ( $\mu$ l) | Solution B (ml) |
|---------|-------------------------------|-----------------|
| 2       | 17                            | 1.7             |
| 4       | 33                            | 3.3             |
| 6       | 50                            | 5.0             |
| 8       | 66                            | 6.6             |
| 10      | 82                            | 8.2             |
| 12      | 100                           | 10.0            |

6. **Wash:** Wash the plate with 1X wash buffer at least 3 times using a wash bottle with manifold or an automated plate washer. Empty the plate by inverting it and blotting on a paper towel to remove excess liquid. *Do not allow the plate to dry out.*
7. **Add Streptavidin Peroxidase Solution:** Prepare streptavidin peroxidase solution with Streptavidin Peroxidase Dilution Buffer (Solution D) as shown in the following table. Add 100  $\mu$ l of streptavidin peroxidase solution to each well and incubate at room temperature for 30 minutes.

| Strip # | Streptavidin Peroxidase ( $\mu$ l) | Solution D (ml) |
|---------|------------------------------------|-----------------|
| 2       | 8                                  | 1.7             |
| 4       | 17                                 | 3.3             |
| 6       | 25                                 | 5.0             |
| 8       | 33                                 | 6.6             |
| 10      | 42                                 | 8.2             |
| 12      | 50                                 | 10.0            |

8. **Wash:** Wash the plate with 1X wash buffer at least 3 times using a wash bottle with manifold or an automated plate washer. Empty the plate by inverting it and blotting on a paper towel to remove excess liquid. *Do not allow the plate to dry out.*
9. **Add TMB Solution:** Use new tubes when preparing TMB solution. Just prior to use, prepare TMB solution with Chromogen Dilution Buffer as shown in the following table. Add 100  $\mu$ l of TMB solution to each well immediately after washing the plate and incubate for 25 minutes at room temperature


| Strip # | TMB ( $\mu$ l) | Chromogen Dilution Buffer (ml) |
|---------|----------------|--------------------------------|
| 2       | 34             | 1.7                            |
| 4       | 66             | 3.3                            |
| 6       | 100            | 5.0                            |
| 8       | 132            | 6.6                            |
| 10      | 164            | 8.2                            |
| 12      | 200            | 10.0                           |

10. **Stop:** Stop the reaction with 50  $\mu$ l of 2N Sulfuric Acid (Stop Solution) to each well.
11. **Read Plate:** Read the OD values at 450 nm. If the OD values of samples are greater than the OD values of the highest standard, re-assay the samples at a higher dilution. A 630 nm filter can be used as a reference.

## CALCULATING RESULTS

1. Average the duplicate OD values for the standards, blanks (B), and test samples.
2. Subtract the averaged blank OD values from the averaged OD values of the standards and test samples.
3. Plot the OD values of the standards against the ng/ml of antibody standard. Using a log/log plot will linearize the data. Figure 1 shows an example of a standard curve where the standard range is 0.8 to 50 ng/ml.
4. The ng/ml of antigen in test samples can be calculated using regression analysis. Multiply it by the sample dilution factor to obtain the mMCP-1 concentration (ng/ml) in the original test samples.

Figure 1 - A Typical Standard Curve for the mouse MCPT-1 Detection ELISA Kit



## VALIDATION DATA

Table 1 - Reproducibility Data for the mouse MCPT-1 Detection ELISA Kit

| Test               | 1.6 ng/ml | 6.3 ng/ml | 25 ng/ml |
|--------------------|-----------|-----------|----------|
| Intra-Assay CV (%) | 6.7       | 4.4       | 3.1      |
| Inter-Assay CV (%) | 2.6       | 1.7       | 8.0      |
| Spike Test* (%)    | 91%       | 108%      | 99%      |

\*Known amounts of mouse MCP-1 were added to samples and then diluted with Sample/Standard/Secondary Antibody Dilution Buffer to assay MCPT-1 by ELISA.

## TROUBLESHOOTING

For frequently asked questions about assays and ELISAs, please see Chondrex, Inc.'s [ELISA FAQ](#) for more information.

---

## REFERENCES

1. S. J. Galli, M. Tsai, IgE and mast cells in allergic disease. *Nat. Med* **18**, 693–704 (2012).
2. G.P. Pejler, E. Rönnberg, I. Waern, S. Wernersson, Mast cell proteases: multifaceted regulators of inflammatory disease. *Blood* **115** (24): 4981–4990 (2010).
3. H. Trong, G. Newlands, H. Miller, H. Charbonneau, H. Neurath, R. Woodbury, Amino acid sequence of a mouse mucosal mast cell protease. *Biochemistry* **28**(1):391-5 (1989).
4. J. Wastling, C. Scudamore, E. Thornton, G. Newlands, H. Miller, Constitutive expression of mouse mast cell protease-1 in normal BALB/c mice and its up-regulation during intestinal nematode infection. *Immunology*, **90**(2), 3 (1997).
5. P. Knight, S. Wright, C. Lawrence, Y. Paterson, H. Miller, Delayed expulsion of the nematode *Trichinella spiralis* in mice lacking the mucosal mast cell-specific granule chymase, mouse mast cell protease-1. *J Exp Med*, **192**(12), (2000).
6. K. Vaali, T. Puumalainen, M. Lehto, H. Wolff, H. Rita, H. Alenius, T. Palosuo, Murine model of food allergy after epicutaneous sensitization: role of mucosal mast cell protease-1. *Scandinavian journal of gastroenterology*, **41**(12), 1405-1 (2006).
7. J. McDermott, R. Bartram, P. Knight, H. Miller, D. Garrod, R. Grecis, Mast cells disrupt epithelial barrier function during enteric nematode infection. *Proceedings of the National Academy of Sciences*, **100**(13), 7761-7766 (2003).