

## Using ChonBlock™ for Anti-Antigen Antibody ELISAs

*For Research Use Only - Not Human or Therapeutic Use*

1. **Coat Plates with Antigen:** Dissolve antigen at 5 - 10 µg/ml in PBS with 0.05% Azide. Add 100 µl of the antigen solution to ELISA plates and incubate at 4°C overnight. Wash the antigen-coated plates with endotoxin-free distilled water 3 times using a wash bottle with manifold. Empty the plate by inverting it and blotting on a paper towel to remove excess liquid. *Do not allow the plate to dry out.*
2. **Add Blocking Buffer:** Add 100 µl of the ChonBlock™ Blocking/Sample Dilution Buffer to each well and incubate for 1 hour at room temperature.
3. **Prepare Sample Dilutions:** Dilute 10 µl of serum sample with 990 µl of ChonBlock™ Blocking/Sample Dilution Buffer (1:100 dilution) and keep it as a stock solution for future assays. Then, further dilute the sample with ChonBlock™ Blocking/Sample Dilution Buffer depending on the antibody levels. For example, take 200 µl of the sample stock solution and mix with 200 µl of ChonBlock™ Blocking/Sample Dilution Buffer to make a 1:200 dilution.

NOTE: Chondrex, Inc. recommends running a preliminary assay using various dilutions of sera (1:200, 1:1,000, 1:5,000) in order to determine the optimal dilution of the samples, especially before assaying a large number of samples.

4. **Wash:** Wash the plate with a wash buffer (PBS with 0.05% Tween 20) at least 3 times using a wash bottle with manifold or an automated plate washer. Empty the plate by inverting it and blotting on a paper towel to remove excess liquid. *Do not allow the plate to dry out.*
  5. **Add Samples:** Add 100 µl of ChonBlock™ Blocking/Sample Dilution Buffer (blank) and diluted samples to wells in duplicate according to the desired layout. Incubate at room temperature for 2 hours.
  6. **Wash:** Wash the plate with the wash buffer at least 3 times using a wash bottle with manifold or an automated plate washer. Empty the plate by inverting it and blotting on a paper towel to remove excess liquid. *Do not allow the plate to dry out.*
  7. **Add Detection Antibody:** Dilute appropriate peroxidase-conjugated detection antibodies in 10 ml of ChonBlock™ Detection Antibody Dilution Buffer. Add 100 µl of the detection antibody solution to each well and incubate at room temperature for 1 hour.
- NOTE: Chondrex, Inc. recommends running preliminary assays using a positive control to optimize the detection antibody concentration.
8. **Wash:** Wash the plate with the wash buffer at least 3 times using a wash bottle with manifold or an automated plate washer. Empty the plate by inverting it and blotting on a paper towel to remove excess liquid. *Do not allow the plate to dry out.*
  9. **Add TMB:** Prepare TMB solution just prior to use. Add 100 µl of the TMB solution to all wells immediately after washing the plate and incubate at room temperature for 10 - 30 minutes.
  10. **Stop:** Add 50 ml of 2N sulfuric acid (Stop Solution) to each well.
  11. **Read Plate:** Read the OD values at 450 nm (a 630 nm filter can be used as a reference).