

## Human Anti-Bacteria & Toxins Antibody ELISA Kits

Catalog # 6113-6128

*For Research Use Only - Not Human or Therapeutic Use*

### PRODUCT SPECIFICATIONS

DESCRIPTION:	ELISA kit to quantify human anti-bacteria and toxin IgG and IgA antibodies
FORMAT:	Pre-coated 96-well ELISA Plate with removeable strips
ASSAY TYPE:	Indirect ELISA
ASSAY TIME:	4.5 hours
STANDARD RANGE:	32 units/ml to 0.5 units/ml
NUMBER OF SAMPLES:	Up to 40 (duplicate) diluted samples/plate and up to 20 (duplicate) low dilution samples/plate
SAMPLE TYPES:	Serum and Plasma
RECOMMENDED SAMPLE DILUTIONS:	1:1000 (at least)
CHROMOGEN:	TMB (read at 450 nm)
STORAGE:	-20°C
VALIDATION DATA:	6113: Intra-Assay (3.2-7.2%)/Inter-Assay (3.7-6.5%)/Spiking Test (95-109%) 6114: Intra-Assay (3.3-8.8%)/Inter-Assay (4.2-9.8%)/Spiking Test (94-103%) 6115: Intra-Assay (1.2-4.7%)/Inter-Assay (6-9.8%)/Spiking Test (90-107%) 6116: Intra-Assay (1.1-5.3%)/Inter-Assay (2-10%)/Spiking Test (88-94%) 6117: Intra-Assay (2.5-8.8%)/Inter-Assay (5.8-9%)/Spiking Test (94-96%) 6118: Intra-Assay (1.1-8.9%)/Inter-Assay (4.5-9.9%)/Spiking Test (105-108%) 6119: Intra-Assay (1.9-7.3%)/Inter-Assay (6-9%)/Spiking Test (95-110%) 6120: Intra-Assay (3.3-9%)/Inter-Assay (1.5-10.9%)/Spiking Test (89-91%) 6121: Intra-Assay (1.3-7.3%)/Inter-Assay (4.4-9.2%)/Spiking Test (92-98%) 6122: Intra-Assay (0.6-6%)/Inter-Assay (6-9.4%)/Spiking Test (93-102%) 6123: Intra-Assay (1.2-8.9%)/Inter-Assay (1.7-8.9%)/Spiking Test (95-108%) 6124: Intra-Assay (3.5-7.3%)/Inter-Assay (4.3-9.7%)/Spiking Test (94-96%) 6125: Intra-Assay (4.2-8%)/Inter-Assay (4-7.8%)/Spiking Test (90-98%) 6126: Intra-Assay (2.9-4.9%)/Inter-Assay (4-8.4%)/Spiking Test (91-93%) 6127: Intra-Assay (1.9-7.5%)/Inter-Assay (6.5-8.3%)/Spiking Test (103-110%) 6128: Intra-Assay (5.8-9%)/Inter-Assay (3.4-9.8%)/Spiking Test (93-102%)
NOTES:	N/A

## Human Anti-Bacteria & Toxins Antibody ELISA Kits

Catalog # 6113-6128

*For Research Use Only - Not Human or Therapeutic Use*

### INTRODUCTION

A growing body of research indicates an association between intestinal bacteria and autoimmune diseases. More specifically, dysbiosis or imbalance of intestinal bacterial flora may contribute to the pathogenesis of Rheumatoid Arthritis (RA), as indicated by many studies on intestinal microbes (1-5). Furthermore, recent studies have also suggested a possible link between RA and periodontal diseases caused by *Porphyromonas gingivalis* (*P. gingivalis*) and *Aggregatibacter actinomycetemcomitans* (*A. actinomycetemcomitans*) (6-9). Normally, intestinal bacteria do not affect the host's health, but under certain conditions, may overcome the host's defenses and exert pathogenic effects. Examples of these include immunosenescence, gastrointestinal disorders such as constipation and diarrhea, or other events such as physical and psychological stress (10-17). It is therefore important to consider these risk factors when studying the etiology of autoimmune diseases. Intestinal bacterial imbalance may increase the levels of pathogenic substances in the gastrointestinal lumen and high mucosal permeability may increase translocation of potential pathogenic agents into the circulatory system. Consequently, pathogens can overwhelm the host's defense functions and cause chronic health problems that can evolve into autoimmune disorders (18, 19).


To facilitate and promote studies that determine immune responses to environmental agents in humans, Chondrex, Inc. provides ELISA kits for assaying human serum antibodies against a variety of potential pathogenic and non-pathogenic environmental agents, of which all humans may be universally exposed to during their lifetime. For more information, please contact us at [support@chondrex.com](mailto:support@chondrex.com). These ELISA kits employ ChonBlock™ (Cat # 9068 and Cat # 90681) assay buffers. ChonBlock™ eliminates non-specific reactions involved in the indirect ELISA, especially false positive reactions caused by hydrophobic binding of immunoglobulins in sample specimens to ELISA plates as reported in detail (20-22).

### PLATE COATING AND SETUP

Antigens	IgG	IgA
<i>E. coli</i> (O111:B4) Lipopolysaccharide (LPS)	Red - 6113	Red - 6114
<i>E. coli</i> (O111:B4)	Yellow - 6115	Yellow - 6116
<i>P. gingivalis</i> LPS	Pink - 6117	Pink - 6118
<i>P. gingivalis</i>	Orange - 6119	Orange - 6120
<i>Lactobacillus casei</i>	Purple - 6121	Purple - 6122
Peptidoglycan Polysaccharide (PG-PS)	Blue - 6123	Blue - 6124
Salmonella	Green - 6125	Green - 6126
Yeast Extract	Gray - 6127	Gray - 6128



## ASSAY OUTLINE



## NOTES BEFORE USING ASSAY

NOTE 1: It is recommended that the standard and samples be run in duplicate.

NOTE 2: Warm up all buffers to room temperature before use.

NOTE 3: Crystals may form in Wash Buffer, 20X when stored at cold temperatures. If crystals have formed, warm the wash buffer by placing the bottle in warm water until crystals are completely dissolved.

NOTE 4: Measure exact volume of buffers using a serological pipet, as extra buffer is provided.

NOTE 5: Cover the plate with plastic wrap or a plate sealer after each step to prevent evaporation from the outside wells of the plate.

NOTE 6: For partial reagent use, please see the assay protocol's corresponding step for the appropriate dilution ratio. For example, if the protocol dilutes 50 µl of a stock solution in 10 ml of buffer for 12 strips, then for 6 strips, dilute 25 µl of the stock solution in 5 ml of buffer. Partially used stock reagents may be kept in their original vials and stored at -20°C for use in a future assay.


NOTE 7: This kit contains animal components from non-infectious animals and should be treated as potential biohazards in use and for disposal.

## ASSAY PROCEDURE

- Add Blocking Buffer:** Add 100  $\mu$ l of the Blocking Buffer (Solution A) to each well and incubate for 1 hour at room temperature.

NOTE: If a sample with a dilution of 1:1,000 or less is assayed, antigen non-coated strips should be used. Solution A must be added to the non-coated wells without prior washing because any contaminants in the vessel containing the washing buffer will bind to the strips. For example, add 100  $\mu$ l of Solution A to the antigen-coated strips (S1) and the corresponding uncoated strips (SB1). Incubate for 1 hour at room temperature.

- Prepare Standard Dilutions:** Dissolve one vial of standard (32 units/vial) with 1 ml Standard/Sample dilution buffer (Solution B) to make a 32 units/ml stock standard solution. Prepare standard serial dilutions by mixing 250  $\mu$ l of the stock standard solution with 250  $\mu$ l of Solution B (16 units/ml). Repeat this procedure to make 8, 4, 2, 1, and 0.5 units/ml standard solutions for a total of 7 serial standard dilutions. Keep the remaining 32 units/ml stock standard solution at -20°C for future assays.



- Prepare Sample Dilutions:** Add 10  $\mu$ l of serum sample to 990  $\mu$ l of Solution B (1:100) and keep it as a stock solution for future assays. Then, further dilute the sample with Solution B depending on the antibody levels. For example, take 50  $\mu$ l of the sample stock solution and mix with 450  $\mu$ l of solution B to make a 1:1,000 dilution. If it is necessary to assay antibodies at a low dilution (less than 1:1,000) due to low antibody levels, antigen uncoated control strips will be necessary. Please contact [support@chondrex.com](mailto:support@chondrex.com) for more information.

NOTE 1: Chondrex, Inc. recommends running a preliminary assay using various dilutions of sera (1:1,000, 1:4,000, 1:16,000) in order to determine the optimal dilution of your samples, especially before assaying a large number of samples.

NOTE 2: The following table shows reference antibody levels from 136 Japanese patients with different diseases. The antibody levels will depend on patient groups.

Antigens	IgG (x1000 units/ml)		IgA (x1000 units/ml)	
	Average	STEDV	Average	STEDV
<i>E.coli</i>	38.7	58.1	15.9	37.9
<i>P.gingivalis</i>	12.4	22.7	4.7	4.9
<i>Lactobacillus casei</i>	30.4	62.2	5.3	11.0
Salmonella	43.0	49.3	19.4	43.8
Yeast Extract	51.2	90.3	6.4	9.9
<i>E.coli</i> -LPS	2.6	3.4	2.1	2.8
PG-LPS	48.1	71.5	2.0	2.1
PG-PS	160.0	154.0	16.6	18.8

4. **Dilute Wash Buffer:** Dilute 50 ml of 20X wash buffer in 950 ml of distilled water (1X wash buffer). Wash the plate with 1X wash buffer at least 3 times using a wash bottle with manifold or an automated plate washer. Empty the plate by inverting it and blotting on a paper towel to remove excess liquid. *Do not allow the plate to dry out.*
5. **Add Standards and Samples:** Add 100  $\mu$ l of standards, Solution B (blank), and samples to wells in duplicate according to the desired layout. Incubate at room temperature for 2 hours.

NOTE: If a sample with a dilution of 1:1,000 or less is assayed, add 100  $\mu$ l of the diluted samples to the antigen-coated strips (S1) and the corresponding uncoated strips (SB1).

6. **Wash:** Wash the plate with 1X wash buffer at least 3 times using a wash bottle with manifold or an automated plate washer. Empty the plate by inverting it and blotting on a paper towel to remove excess liquid. Do not allow the plate to dry out.
7. **Add Secondary Antibody:** Dilute one vial of Secondary Antibody with 10 ml of Secondary Antibody Dilution Buffer (Solution C). Add 100  $\mu$ l of secondary antibody solution to each well and incubate at room temperature for 1 hour.

Strip #	2 <sup>nd</sup> Antibody ( $\mu$ l)	Solution C (ml)
2	8	1.7
4	17	3.3
6	25	5.0
8	33	6.6
10	42	8.2
12	50	10.0

8. **Wash:** Wash the plate with 1X wash buffer at least 3 times using a wash bottle with manifold or an automated plate washer. Empty the plate by inverting it and blotting on a paper towel to remove excess liquid. Do not allow the plate to dry out.
9. **Add TMB:** Use new tubes when preparing TMB. Dilute one vial of TMB with 10 ml of Chromogen Dilution Buffer just prior to use. Add 100  $\mu$ l of the TMB solution to all wells immediately after washing the plate and incubate at room temperature for 25 minutes.

Strip #	TMB ( $\mu$ l)	Chromogen Dilution Buffer (ml)
2	34	1.7
4	66	3.3
6	100	5.0
8	132	6.6
10	164	8.2
12	200	10.0

10. **Stop:** Add 50  $\mu$ l of 2N sulfuric acid (Stop Solution) to each well.
11. **Read Plate:** Read the OD values at 450 nm (a 630 nm filter can be used as a reference) within 5 minutes after added Stop Solution. If the OD values of the samples are greater than the OD values of the highest standard, re-assay the samples at a higher dilution.


## CALCULATING RESULTS

1. Average the duplicate OD values for the standards, blanks (B) and test samples in uncoated wells and coated wells.
2. Subtract the blank (B) values from the averaged OD values of the standards and test samples in uncoated wells and coated wells.

NOTE: Individual antigens have unique background values. Therefore, blank wells should be used for each different antigen.

3. Subtract the OD values of samples tested in uncoated wells (background values) from their counterpart OD values in coated wells from Step 2 to eliminate values associated with non-specific reactions.
4. Plot the OD values of standards against the units/ml of antibody standard. Using a log/log plot will linearize the data. Figure 1 shows a representative experiment where the standard range is from 0.5 to 32 units/ml. The units/ml of antibody in test samples can be calculated using regression analysis.

Figure 1 - A Typical Standard Curve for a Human-Anti-Bacteria and Toxin Antibody ELISA Kit



## VALIDATION DATA

Table 1 - Reproducibility data for the Anti-*E. coli* (O111:B4) LPS Antibody ELISA Kit

Test At	IgG (6113)			IgA (6114)		
	1 unit/ml	4 units/ml	16 units/ml	1 unit/ml	4 units/ml	16 units/ml
Intra-Assay CV (%)	3.2	6.2	7.2	8.8	5.0	3.3
Inter-Assay CV (%)	3.7	6.3	6.5	9.8	9.3	4.2
Spiking Test*	97%	109%	95%	103%	97%	94%

Table 2 - Reproducibility data for the Anti-*E. coli* (O111:B4) Antibody ELISA Kit

Test At	IgG (6115)			IgA (6116)		
	1 unit/ml	4 units/ml	16 units/ml	1 unit/ml	4 units/ml	16 units/ml
Intra-Assay CV (%)	4.0	4.7	1.2	1.1	3.8	5.3
Inter-Assay CV (%)	9.8	8.0	6.0	3.8	2.0	10.0
Spiking Test*	95%	107%	90%	88%	94%	91%



Table 3 - Reproducibility data for the Anti-*P. gingivalis* LPS Antibody ELISA Kit

Test At	IgG (6117)			IgA (6118)		
	1 unit/ml	4 units/ml	16 units/ml	1 unit/ml	4 units/ml	16 units/ml
Intra-Assay CV (%)	2.5	8.8	2.8	8.9	1.1	1.9
Inter-Assay CV (%)	9.0	5.8	7.7	9.9	4.5	6.3
Spiking Test*	94%	96%	95%	105%	108%	105%

Table 4 - Reproducibility data for the Anti-*P. gingivalis* Antibody ELISA Kit

Test At	IgG (6119)			IgA (6120)		
	1 unit/ml	4 units/ml	16 units/ml	1 unit/ml	4 units/ml	16 units/ml
Intra-Assay CV (%)	7.3	1.9	2.8	3.3	5.3	9.0
Inter-Assay CV (%)	6.0	6.3	9.0	10.9	1.5	5.8
Spiking Test*	95%	95%	110%	91%	89%	91%

Table 5 - Reproducibility data for the Anti-*Lactobacillus casei* Antibody ELISA Kit

Test At	IgG (6121)			IgA (6122)		
	1 unit/ml	4 units/ml	16 units/ml	1 unit/ml	4 units/ml	16 units/ml
Intra-Assay CV (%)	3.9	7.3	1.3	6.0	2.6	0.6
Inter-Assay CV (%)	9.2	4.4	7.3	6.0	9.2	9.4
Spiking Test*	92%	94%	98%	96%	102%	93%

Table 6 - Reproducibility data for the Anti-PG-PS Antibody ELISA Kit

Test At	IgG (6123)			IgA (6124)		
	1 unit/ml	4 units/ml	16 units/ml	1 unit/ml	4 units/ml	16 units/ml
Intra-Assay CV (%)	6.8	1.2	8.9	7.3	3.5	6.5
Inter-Assay CV (%)	6.1	1.7	8.9	7.5	9.7	4.3
Spiking Test*	108%	101%	95%	94%	96%	96%

Table 7 - Reproducibility data for the Anti-Salmonella Antibody ELISA Kit

Test At	IgG (6125)			IgA (6126)		
	1 unit/ml	4 units/ml	16 units/ml	1 unit/ml	4 units/ml	16 units/ml
Intra-Assay CV (%)	4.2	5.0	8.0	4.9	3.6	2.9
Inter-Assay CV (%)	4.6	4.0	7.8	8.4	4.8	4.0
Spiking Test*	90%	91%	98%	92%	91%	93%

Table 8 - Reproducibility data for the Anti-Yeast Extract Antibody ELISA Kit

Test At	IgG (6127)			IgA (6128)		
	1 unit/ml	4 units/ml	16 units/ml	1 unit/ml	4 units/ml	16 units/ml
Intra-Assay CV (%)	7.5	6.1	1.9	6.6	9.0	5.8
Inter-Assay CV (%)	6.5	6.9	8.3	9.7	3.4	9.8
Spiking Test*	110%	103%	104%	95%	102%	93%

\* Standards were added with known amounts of antibody and then diluted with Sample/Standard Dilution Buffer to assay antibodies by ELISA.

## TROUBLESHOOTING

For frequently asked questions about assays and ELISAs, please see Chondrex, Inc.'s [ELISA FAQ](#) for more information.



## REFERENCES

1. J. Kjeldsen-Kragh, Rheumatoid Arthritis Treated With Vegetarian Diets. *Am J Clin Nutr* **70**, 594S-600S (1999).
2. R. Peltonen, M. Nenonen, T. Helve, O. Hänninen, P. Toivanen, E. Eerola, *et al.*, Faecal Microbial Flora and Disease Activity in Rheumatoid Arthritis During a Vegan Diet. *Br J Rheumatol* **36**, 64-8 (1997).
3. J. Vaahtovuori, E. Munukka, M. Korkeamäki, R. Luukkainen, P. Toivanen, Fecal Microbiota in Early Rheumatoid Arthritis. *J Rheumatol* **35**, 1500-5 (2008).
4. J. Scher, A. Sczesnak, R. Longman, N. Segata, C. Ubeda, *et al.*, Expansion of Intestinal *Prevotella Copri* Correlates With Enhanced Susceptibility to Arthritis. *Elife* **2**, e01202 (2013).
5. Y. Kang, Y. Cai, X. Zhang, X. Kong, J. Su, Altered Gut Microbiota in RA: Implications for Treatment. *Z Rheumatol* **76**, 451-457 (2017).
6. P. de, T. Dietrich, T. McAlindon, Association of Periodontal Disease and Tooth Loss With Rheumatoid Arthritis in the US Population. *J Rheumatol* **35**, 70-6 (2008).
7. F. Liao, Z. Li, Y. Wang, B. Shi, Z. Gong, X. Cheng, *et al.*, *Porphyromonas Gingivalis* May Play an Important Role in the Pathogenesis of Periodontitis-Associated Rheumatoid Arthritis. *Med Hypotheses* **72**, 732-5 (2009).
8. Y. Chou, K. Lai, D. Chen, C. Lin, H. Chen, Rheumatoid Arthritis Risk Associated With Periodontitis Exposure: A Nationwide, Population-Based Cohort Study. *PLoS One* **10**, e0139693 (2015).
9. A. Yoshida, Y. Nakano, Y. Yamashita, T. Oho, H. Ito, *et al.*, Immunodominant Region of *Actinobacillus Actinomycetemcomitans* 40-kilodalton Heat Shock Protein in Patients With Rheumatoid Arthritis. *J Dent Res* **80**, 346-50 (2001).
10. D. Schmucker, K. Thoreux, R. Owen, Aging Impairs Intestinal Immunity. *Mech Ageing Dev* **122**, 1397-411 (2001).
11. K. Tiihonen, A. Ouwehand, N. Rautonen, Human Intestinal Microbiota and Healthy Ageing. *Ageing Res Rev* **9**, 107-16 (2010).
12. A. Larbi, P. Rymkiewicz, A. Vasudev, I. Low, N. Shadan, S. Mustafah, *et al.*, The Immune System in the Elderly: A Fair Fight Against Diseases? *Ageing Health* **9**, 35-47 (2013).
13. L. Ferrucci, A. Corsi, F. Lauretani, S. Bandinelli, B. Bartali, *et al.*, The Origins of Age-Related Proinflammatory State. *Blood* **105**, 2294-9 (2005).
14. Khalif, E. Quigley, E. Konovitch, I. Maximova, Alterations in the Colonic Flora and Intestinal Permeability and Evidence of Immune Activation in Chronic Constipation. *Dig Liver Dis* **37**, 838-49 (2005).
15. K. Katayama, T. Matsuno, T. Waritani, K. Terato, H. Shionoya, Supplemental Treatment of Rheumatoid Arthritis With Natural Milk Antibodies Against Enteromicrobes and Their Toxins: Results of an Open-Labelled Pilot Study. *Nutr J* **10**, 2 (2011).
16. P. Anderlik, I. Szeri, Z. Bános, Z. Barna, Bacterial Translocation After Cold Stress in Young and Old Mice. *Acta Microbiol Hung* **37**, 289-94 (1990).
17. G. Selkirk, T. McLellan, H. Wright, S. Rhind, Mild Endotoxemia, NF-kappaB Translocation, and Cytokine Increase During Exertional Heat Stress in Trained and Untrained Individuals. *Am J Physiol Regul Integr Comp Physiol* **295**, R611-23 (2008).
18. K. Terato, C. Do, H. Shionoya, Slipping Through the Cracks: Linking Low Immune Function and Intestinal Bacterial Imbalance to the Etiology of Rheumatoid Arthritis. *Autoimmune Dis* **2015**, 636207 (2015).
19. K. Terato, T. Waritani, R. Fukai, H. Shionoya, H. Itoh, K. Katayama, *et al.*, Contribution of Bacterial Pathogens to Evoking Serological Disease Markers and Aggravating Disease Activity in Rheumatoid Arthritis. *PLoS One* **13**, e0190588 (2018).
20. K. Terato, C. Do, D. Cutler, T. Waritani, H. Shionoya, Preventing Intense False Positive and Negative Reactions Attributed to the Principle of ELISA to Re-Investigate Antibody Studies in Autoimmune Diseases. *J Immunol Methods* **407**, 15-25 (2014).
21. K. Terato, C. Do, J. Chang, T. Waritani, Preventing Further Misuse of the ELISA Technique and Misinterpretation of Serological Antibody Assay Data. *Vaccine* **34**, 4643-4644 (2016).
22. T. Waritani, J. Chang, B. McKinney, K. Terato, An ELISA Protocol to Improve the Accuracy and Reliability of Serological Antibody Assays. *MethodsX* **4**, 153-165 (2017).