

Mouse anti-*Escherichia coli* (O111:B4) IgM and IgG3 Antibody ELISA Kits

Catalog # 6209 and 6212

For Research Use Only - Not Human or Therapeutic Use

PRODUCT SPECIFICATIONS

DESCRIPTION:	ELISA kits to quantify mouse anti- <i>E. coli</i> IgM and IgG3 antibodies 6209: Mouse Anti- <i>E. coli</i> IgM Antibody ELISA Kit 6212: Mouse Anti- <i>E. coli</i> IgG3 Antibody ELISA Kit
FORMAT:	Pre-coated 96-well ELISA Plate with removeable strips
ASSAY TYPE:	Indirect ELISA
ASSAY TIME:	4.5 hours
STANDARD RANGE:	500 ng/ml to 8 ng/ml
NUMBER OF SAMPLES:	Up to 40 (duplicate) diluted samples/plate and up to 20 (duplicate) low dilution samples/plate
SAMPLE TYPES:	Serum and Plasma
RECOMMENDED SAMPLE DILUTIONS:	1:200 (at least)
CHROMOGEN:	TMB (read at 450 nm)
STORAGE:	-20°C
VALIDATION DATA:	6209: Intra-Assay (6-7.6%)/Inter-assay (3.9-10%)/Spiking Test (104-108%) 6212: Intra-Assay (2.8-4.9%)/Inter-assay (5.7-9.7%)/Spiking Test (94-108%)
NOTES:	N/A

Mouse anti-*Escherichia coli* (O111:B4) IgM and IgG3 Antibody ELISA Kits

Catalog # 6209 and 6212

For Research Use Only - Not Human or Therapeutic Use

INTRODUCTION

Recent studies indicate environmental factors, especially intestinal microbiota and their toxins, may play pathogenic roles in autoimmune diseases such as rheumatoid arthritis (RA) (1-7), inflammatory bowel diseases (IBDs) (8, 9), systemic lupus erythematosus (SLE) (10), and other chronic disorders (11-13). Increased translocation of bacteria and bacterial toxins associated with high mucosal permeability and low immune function may play roles in the primary and common pathogenesis of these autoimmune disorders (14). Mice are ideal experimental animals to study this hypothesis due to the variety in genetic backgrounds, strains, and housing conditions (germ free, gnotobiotic, specific pathogen free, or conventional), which can affect susceptibility to potential pathogens. The variability of susceptibility to pathogens may contribute to the maturity of the immune system, depending on the environment.

To facilitate and promote studies that determine immune responses to environmental agents, Chondrex, Inc. provides mouse antibody ELISA kits against *Escherichia coli* O111:B4 (*E. coli*). In addition, Chondrex, Inc. also provides mouse antibody ELISA kits against lipopolysaccharide from *E. coli* O111:B4, *Staphylococcus aureus*, and Staphylococcal enterotoxins A & B. For further requests or consultation, please contact support@chondrex.com or visit www.chondrex.com.

Catalog

Mouse Anti-*E. coli* IgM Antibody ELISA Kit (Cat # 6209)
Mouse Anti-*E. coli* IgG3 Antibody ELISA Kit (Cat # 6212)

KIT COMPONENTS

Item	Quantity	Amount	Storage
Standard Antibody IgM (62096) IgG3 (62126)	1 vial	500 ng, lyophilized	-20°C
Secondary Antibody (peroxidase-conjugated polyclonal antibodies) IgM (62092) IgG3 (62122)	2 vials	50 µl/vial	-20°C
Solution A - Blocking Buffer (61026)	1 bottle	10 ml	-20°C
Solution B - Sample/Standard Dilution Buffer (61027)	1 bottle	50 ml	-20°C
Solution C - Secondary Antibody Dilution Buffer (61025)	1 bottle	20 ml	-20°C
TMB Solution (90023)	2 vials	0.2 ml/vial	-20°C
Chromogen Dilution Buffer (90022)	1 bottle	20 ml	-20°C
Stop Solution - 2N Sulfuric Acid (9016)	1 bottle	10 ml	-20°C
Wash Buffer, 20X (9005)	1 bottle	50 ml	-20°C
<i>E. coli</i> O111:B4-coated 8-Well Strips (Yellow)	12 each	8-well strips	-20°C

An antigen un-coated plate (Cat #9026) is not included. Please contact support@chondrex.com to place an order.

ASSAY OUTLINE






PLATE MAPPING

Map the plate based on the number of samples and sample dilution. For example, if sample dilution is more than 1:200, it is not necessary to run antigen uncoated wells, but if sample dilution is less than 1:200, it is necessary to run antigen uncoated wells to determine the background noise (BG) reaction OD values of individual samples.

Standard Layout of Antigen-coated Plate



Standard Layout of Antigen-coated and Uncoated Plate



NOTES BEFORE USING ASSAY

NOTE 1: It is recommended that the standard and samples be run in duplicate.

NOTE 2: Warm up all buffers to room temperature before use.

NOTE 3: Crystals may form in Wash Buffer, 20X when stored at cold temperatures. If crystals have formed, warm the wash buffer by placing the bottle in warm water until crystals are completely dissolved.

NOTE 4: Measure exact volume of buffers using a serological pipet, as extra buffer is provided.

NOTE 5: Cover the plate with plastic wrap or a plate sealer after each step to prevent evaporation from the outside wells of the plate.

NOTE 6: For partial reagent use, please see the assay protocol's corresponding step for the appropriate dilution ratio. For example, if the protocol dilutes 50 μ l of a stock solution in 10 ml of buffer for 12 strips, then for 6 strips, dilute 25 μ l of the stock solution in 5 ml of buffer. Partially used stock reagents may be kept in their original vials and stored at -20°C for use in a future assay.


NOTE 7: This kit contains animal components from non-infectious animals and should be treated as potential biohazards in use and for disposal.

ASSAY PROCEDURE

- Add Blocking Buffer:** Add 100 μ l of the Blocking Buffer (Solution A) to each well and incubate for 1 hour at room temperature.

NOTE: If a sample with a dilution of 1:200 or less is assayed, antigen non-coated strips should be used. Solution A must be added to the non-coated wells without prior washing because any contaminants in the vessel containing the washing buffer will bind to the strips. For example, add 100 μ l of Solution A to the antigen-coated strips (S1) and the corresponding uncoated strips (SB1). Incubate for 1 hour at room temperature.

- Prepare Standard Dilutions:** Dissolve one vial of standard (500 ng/vial) with 1 ml Standard/Sample dilution buffer (Solution B) to make a 500 ng/ml stock standard solution. Prepare standard serial dilutions by mixing 250 μ l of the stock standard solution with 250 μ l of Solution B (250 ng/ml). Repeat this procedure to make 125, 63, 31, 16, and 8 ng/ml standard solutions for a total of 7 serial standard dilutions. Keep the remaining 500 ng/ml stock standard solution at -20°C for future assays.



- Prepare Sample Dilutions:** Add 10 μ l of serum sample to 990 μ l of Solution B (1:100) and keep it as a stock solution for future assays. Then, further dilute the sample with Solution B depending on the antibody levels. For example, take 200 μ l of the sample stock solution and mix with 200 μ l of solution B to make a 1:200 dilution. If it is necessary to assay antibodies at a low dilution (less than 1:200) due to low antibody levels, antigen uncoated control strips will be necessary. Please contact support@chondrex.com for more information.

NOTE: Chondrex, Inc. recommends running a preliminary assay using various dilutions of sera (1:200, 1:1,000, 1:5,000) in order to determine the optimal dilution of your samples, especially before assaying a large number of samples.

- Dilute Wash Buffer:** Dilute 50 ml of 20X wash buffer in 950 ml of distilled water (1X wash buffer). Wash the plate with 1X wash buffer at least 3 times using a wash bottle with manifold or an automated plate washer. Empty the plate by inverting it and blot on a paper towel to remove excess liquid. Do not allow the plate to dry out.
- Add Standards and Samples:** Add 100 μ l of standards, Solution B (blank), and samples to wells in duplicate. Incubate at room temperature for 2 hours.

NOTE: If a sample with a dilution of 1:200 or less is assayed, add 100 μ l of the diluted samples to the antigen-coated strips (S1) and the corresponding uncoated strips (SB1).

- Wash:** Wash the plate with 1X wash buffer at least 3 times using a wash bottle with manifold or an automated plate washer. Empty the plate by inverting it and blot on a paper towel to remove excess liquid. Do not allow the plate to dry out.
- Add Secondary Antibody:** Dilute one vial of Secondary Antibody with 10 ml of Secondary Antibody Dilution Buffer (Solution C). Add 100 μ l of secondary antibody solution to each well and incubate at room temperature for 1 hour.

Strip #	2 nd Antibody (μ l)	Solution C (ml)
2	8	1.7
4	17	3.3
6	25	5.0
8	33	6.6
10	42	8.2
12	50	10.0

- Wash:** Wash the plate with 1X wash buffer at least 3 times using a wash bottle with manifold or an automated plate washer. Empty the plate by inverting it and blot on a paper towel to remove excess liquid. Do not allow the plate to dry out.
- Add TMB:** Use new tubes when preparing TMB. Dilute one vial of TMB with 10 ml of Chromogen Dilution Buffer just prior to use. Add 100 μ l of the TMB solution to all wells immediately after washing the plate and incubate at room temperature for 25 minutes.


Strip #	TMB (μ l)	Chromogen Dilution Buffer (ml)
2	34	1.7
4	66	3.3
6	100	5.0
8	132	6.6
10	164	8.2
12	200	10.0

- Stop:** Add 50 μ l of 2N sulfuric acid (Stop Solution) to each well.
- Read Plate:** Read the OD values at 450 nm (a 630 nm filter can be used as a reference) within 5 minutes after added Stop Solution. If the OD values of the samples are greater than the OD values of the highest standard, re-assay the samples at a higher dilution.

CALCULATING RESULTS

1. Average the duplicate OD values for the standards, blanks (B) and test samples in uncoated wells and coated wells.
2. Subtract the blank (B) values from the averaged OD values of the standards and test samples in uncoated wells and coated wells.
3. Subtract the OD values of samples tested in uncoated wells (background values) from their counterpart OD values in coated wells from Step 2 to eliminate values associated with non-specific reactions.
4. Plot the OD values of standards against the ng/ml of antibody standard. Using a log/log plot will linearize the data. Figure 1 shows a representative experiment where the standard range is from 8 to 500 ng/ml. The ng/ml of antibody in test samples can be calculated using regression analysis.

Figure 1 -A Typical Standard Curve for the Mouse Anti-*E. coli* IgG3 Antibody ELISA Kit



VALIDATION DATA

Table 1 - Reproducibility data for the Anti-*E. coli* (O111:B4) IgM Antibody Assay Kit

Test	15.6 ng/ml	62.5 ng/ml	250 ng/ml
Intra-Assay CV (%)	7.6	6.0	6.9
Inter-Assay CV (%)	10.0	8.8	3.9
Spike Test* (%)	108%	104%	107%

Table 2 - Reproducibility data for the Anti-*E. coli* (O111:B4) IgG3 Antibody Assay Kit

Test	15.6 ng/ml	62.5 ng/ml	250 ng/ml
Intra-Assay CV (%)	4.9	4.7	2.8
Inter-Assay CV (%)	5.7	9.7	7.5
Spike Test* (%)	108%	105%	94%

* Known amounts of antibody were added to standard amounts and then diluted with Sample/Standard Dilution Buffer to assay anti-*E. coli* (O111:B4) antibodies by ELISA.

TROUBLESHOOTING

For frequently asked questions about assays and ELISAs, please see Chondrex, Inc.'s [ELISA FAQ](#) for more information.

REFERENCES

1. S. Aoki, K. Yoshikawa, T. Yokoyama, T. Nonogaki, S. Iwasaki, *et al.*, Role of Enteric Bacteria in the Pathogenesis of Rheumatoid Arthritis: Evidence for Antibodies to Enterobacterial Common Antigens in Rheumatoid Sera and Synovial Fluids. *Ann Rheum Dis* **55**, 363-9 (1996).
2. I. van der Heijden, B. Wilbrink, I. Tchetverikov, I. Schrijver, L. Schouls, *et al.*, Presence of Bacterial DNA and Bacterial Peptidoglycans in Joints of Patients With Rheumatoid Arthritis and Other Arthritides. *Arthritis Rheum* **43**, 593-8 (2000).
3. K. Terato, X. Ye, H. Miyahara, M. Cremer, M. Griffiths, Induction by Chronic Autoimmune Arthritis in DBA/1 Mice by Oral Administration of Type II Collagen and Escherichia Coli Lipopolysaccharide. *Br J Rheumatol* **35**, 828-38 (1996).
4. R. Peltonen, M. Nenonen, T. Helve, O. Hänninen, P. Toivanen, E. Eerola, *et al.*, Faecal Microbial Flora and Disease Activity in Rheumatoid Arthritis During a Vegan Diet. *Br J Rheumatol* **36**, 64-8 (1997).
5. P. Toivanen, Normal Intestinal Microbiota in the Aetiopathogenesis of Rheumatoid Arthritis. *Ann Rheum Dis* **62**, 807-11 (2003).
6. J. Vaahтовuо, E. Munukka, M. Korkeamäki, R. Luukkainen, P. Toivanen, Fecal Microbiota in Early Rheumatoid Arthritis. *J Rheumatol* **35**, 1500-5 (2008).
7. C. Edwards, Commensal Gut Bacteria and the Etiopathogenesis of Rheumatoid Arthritis. *J Rheumatol* **35**, 1477-14797 (2008).
8. D. Shi, J. Das, G. Das, Inflammatory Bowel Disease Requires the Interplay Between Innate and Adaptive Immune Signals. *Cell Res* **16**, 70-4 (2006).
9. S. Nell, S. Suerbaum, C. Josenhans, The Impact of the Microbiota on the Pathogenesis of IBD: Lessons From Mouse Infection Models. *Nat Rev Microbiol* **8**, 564-77 (2010).
10. T. Cavallo, N. Granholm, Bacterial Lipopolysaccharide Induces Long-Lasting IgA Deficiency Concurrently With Features of Polyclonal B Cell Activation in Normal and in Lupus-Prone Mice. *Clin Exp Immunol* **84**, 134-8 (1991).
11. W. Penhale, P. Young, The Influence of the Normal Microbial Flora on the Susceptibility of Rats to Experimental Autoimmune Thyroiditis. *Clin Exp Immunol* **72**, 288-92 (1988).
12. M. Murakami, T. Tsubata, R. Shinkura, S. Nisitani, M. Okamoto, *et al.*, Oral Administration of Lipopolysaccharides Activates B-1 Cells in the Peritoneal Cavity and Lamina Propria of the Gut and Induces Autoimmune Symptoms in an Autoantibody Transgenic Mouse. *J Exp Med* **180**, 111-21 (1994).
13. M. Nymark, P. Pussinen, A. Tuomainen, C. Forsblom, P. Groop, *et al.*, Serum Lipopolysaccharide Activity Is Associated With the Progression of Kidney Disease in Finnish Patients With Type 1 Diabetes. *Diabetes Care* **32**, 1689-93 (2009).
14. K. Terato, C. Do, D. Cutler, T. Waritani, H. Shionoya, Preventing Intense False Positive and Negative Reactions Attributed to the Principle of ELISA to Re-Investigate Antibody Studies in Autoimmune Diseases. *J Immunol Methods* **407**, 15-25 (2014).